

Quality Products that Work for You...

Quality comes first at Stone. We've been making handcrafted products for livestock producers since 1932. Over the years we've learned a great deal about quality. You can see it in the meticulous detail and workmanship of our products, and you can feel it in the way they perform.



We know you're only going to choose Stone products if they perform to your complete satisfaction. Above all, we understand that our products must not only make your job safer, easier and more effective, they must also add to your bottom line.

We value quality and performance just as much as you do. The best way to appreciate the true quality of a Stone product is to experience it for yourself. We take pride in making products that work for you.




STONE
MANUFACTURING & SUPPLY COMPANY

Manufacturing & Supply Company
1212 Kansas Avenue
Kansas City, Missouri 64127
(816) 231-4020 • FAX: (816) 241-3336
www.stonemfg.net

BLB SM-5385 9/02 7.5M

Socorex Syringe

Socorex Syringes and Accessories



STONE
MANUFACTURING & SUPPLY COMPANY

Socorex Syringe

Quality, Precision and Safety

Always give the right dosage with this self-refilling syringe with spring-loaded plunger and 3-way valve. Available in six capacities (0.5ml to 20ml) to meet all your vaccination needs.

- Ergonomic shape, perfect balance for fatigue-free handling.
- Precision volume setting with high dose-to-dose reproducibility.
- Unique hassle-free construction – 6 main parts.
- No O-rings to replace, easy to disassemble, clean and maintain.
- Stainless steel plunger works directly against hard machined glass barrel.
- Fully autoclavable at 248° F.



Easy to disassemble, clean and maintain.
Only 6 main parts and no O-rings on the plunger.

**Only 6
main parts**



Socorex Syringe with Vial Holder



Socorex Syringe with Feed Tube

Socorex Twin Syringe with Feed Tubes



SOCOREX SYRINGES

0.5 ml to 20 ml

Capacity	Volume Range	Increment	SOCOREX SYRINGE with Vial Holder Luer Lock Product #	SOCOREX SYRINGE with Feed Tube Kit Luer Lock Product #	SOCOREX SYRINGE with Feed Tube and Vial Holder Luer Lock Product #	SOCOREX TWIN SYRINGE with Feed Tube Kits Luer Lock Product #
0.5 ml	0.05-0.5 ml	0.1 ml	9300	9400	9500	9100
1 ml	0.1-1 ml	.05 ml	9301	9401	9501	9101
2 ml	0.3-2 ml	.1 ml	9302	9402	9502	9102
5 ml	0.5-5 ml	0.5 ml	9305	9405	9505	9105
10 ml	1-10 ml	1 ml	9310	9410	9510	—
20 ml	5-20 ml	1 ml	—	9420	—	—

Socorex Syringe Accessories

Manifold for Socorex Twin Syringe



9570

Allows for the twin syringe to be used for one needle injection.

Product # 9570 – Manifold for Socorex Twin Syringe

Spare Parts for Valves

Top quality set of spares for valves for all six sizes of Socorex syringes. Spare valve kits include: 2 each, stainless steel balls, springs, washers, 1 gasket.



9700

Product # 9700 – 0.5 ml spare valve kit

Product # 9701 – 1 ml spare valve kit

Product # 9702 – 2 ml spare valve kit

Product # 9705 – 5 ml spare valve kit

Product # 9710 – 10 ml spare valve kit

Product # 9920 – Spare valve springs (package 10)

Product # 9921 – Spare valve stainless steel balls (package 10)

Product # 9922 – Spare valve washers (package 10)

NOTE: Call for spare valve parts for 20 ml syringe.

Vial Holder Kit



9517

Kit includes 2-piece slanted vial holder, locking ring, bottle-lock, tube connector and intake needle.

Product # 9517 – Vial Holder Kit for Socorex

Spare Glass Barrels

0.5 ml, 1 ml, 2 ml, 5 ml, 10 ml and 20 ml for Socorex syringes.

Product # 9600 – 0.5 ml glass barrel

Product # 9601 – 1 ml glass barrel

Product # 9602 – 2 ml glass barrel

Product # 9605 – 5 ml glass barrel

Product # 9610 – 10 ml glass barrel

Product # 9620 – 20 ml glass barrel



9605

Feed Tube Kit



9515

Kit includes 1 tube (39.5"), 1 sinker, 1 aspirating cannula (7.5" Feed Needle) and 1 vent cannula (disposable needle). (Does not fit 20 ml syringe.)

Product # 9515 – Feed Tube Kit for Socorex

Holster for Socorex

Durable, adjustable holster features pockets for syringe and vaccine bottle and a pocket for spare needles, marker or other accessories. Allows one-hand operation of syringe when used with Feed Tube Kit.

Product # 9519 – Holster for Socorex



9519

Feed Tube with Bottle Spike



9825

Use this handy bottle spike for inverted or vertical vials. Use with Socorex or other style syringes.

Product # 9825 – Feed Tube with Bottle Spike

Socorex Syringe Accessories

Syringe Cleaning Brushes




- Product # 9250** – Syringe Valve Cleaning Brush
- Product # 9251** – Syringe Brush, 1 ml barrel
- Product # 9252** – Syringe Brush, 2 ml barrel
- Product # 9255** – Syringe Brush, 5 ml barrel
- Product # 9260** – Syringe Brush, 10 ml barrel

Slap-Shot® Flexible Vaccinator

Vaccinate or administer injectable products without total confinement. Inject animals in chutes, alleyways, crates or in the open.

- Minimal needle breakage.
- Attaches to any standard Luer lock syringe and needle.
- Clear flexible tubing stays flexible.
- Empty, rinse and sterilize to re-use.

Product # 9230 – Slap-Shot, 30 inch



Slap-Shot shown with Socorex syringe



Sterimatic

Sterimatic needle guard and sanitizing system. Needle guard protects user from accidental contact with needle, aids in reducing needle breakage and, when used in conjunction with Stericaps, sanitizes the needle at time of injection and at time of needle withdrawal. Available in two sizes. Each Sterimatic comes packaged with 5 Stericaps. Additional Stericaps are available: Product # 9860.



- Product # 9850** – Sterimatic, for subcutaneous (1/2 inch) injections
- Product # 9855** – Sterimatic, for IM (1.5 inch) injections

Sterimatic Socorex Adaptor

Stainless steel Luer lock adaptor for attaching Sterimatic to Socorex syringe.

Product # 9901 – Sterimatic Socorex Adaptor



Stericap Refill Pack

Stericaps are used in conjunction with the Sterimatic system to clean and sanitize needles at time of injection. Sanitized needles can aid in the reduction of abscesses associated with dirty needles and other foreign materials entering injection site.

Product # 9860 – Stericap Refill Pack, package of 5

